



Trust the leader in medical-grade refrigeration for CDC-compliant solutions

## Does your vaccine storage equipment meet CDC guidelines?

Follett's data logger accessory kit is completely compliant with CDC guidelines for digital data logging. The kit can be installed on any Follett refrigerator or freezer as well as on medical-grade models made by others.

## Follett's digital data logger accessory kit meets or exceeds the following CDC guidelines

- Glycol-encased probe
- Calibrated thermometer with an ILAC-accredited Certificate of Traceability and Calibration Testing
- Accuracy of  $\pm 1$  F ( $\pm 0.5$  C)
- Detachable probe to download temperature data without disturbing the glycol-encased probe
- Digital data logger that records and stores temperature information at frequent intervals for continuous temperature monitoring
- Programmable sampling rate with minimum 1 reading every 15 minutes
- Recorder programmed to stop without overwriting data when memory is full
- Current and min/max temperature displays
- Hi/Lo temperature alarms

Data logger features a touchscreen display



## Kit includes

Dickson® data logger | probe | instructions | AC adaptor | Certificate of Traceability and Calibration Testing | Dicksonware™ software USB flash drive | mounting kit | Velcro® | 60 ml bottle | bottle stand | wire ties

Velcro is a registered trademark of Velcro industries B. V.  
Dickson is a registered trademark of Dickson Corporation.  
Dicksonware is a trademark of Dickson Corporation.  
Follett is a registered trademark of Follett LLC, registered in the US.



Contact your Follett representative, call 800.523.9361, or visit [follettice.com/datalogger](http://follettice.com/datalogger) for more information.